

Open Field Test

DESCRIPTION

The Open Field test is available for the measurement of anxious behavior in mice. Animals are placed in the apparatus and are motivated by their natural exploratory behavior. The time spent in the center and edge of the chamber is recorded. The Animal Phenotyping Core houses two of these units.

FEATURES & SPECIFICATIONS

Clear, sturdy apparatus that is easily removed for cleaning

Coupled with Any-Maze behavioral tracking software for objective data collection and analysis

Apparatus Dimensions

- 40 x 40 x 35 cm

Any-Maze behavioral tracking software

- Version 5.2

PHOTOS

